



**TOPSCIENCE 众山**



# STERILEHYA

**Sterile  
Sodium Hyaluronate**



**SHANDONG TOPSCIENCE BIOTECH CO., LTD.**

Address: #98 Lanshan West Road, Lanshan District,  
Rizhao, Shandong, China

Tel: +86-532-80821298

E-mail: [sales@topscience.cn](mailto:sales@topscience.cn)

Website: [www.topscience.cn](http://www.topscience.cn)

**QINGDAO TRONGEN BIOTECH CO., LTD.**

Address: #330 Songling Road, Laoshan District,  
Qingdao, Shandong, China

Tel: +86-532-80821298

E-mail: [sales@trongenbio.com](mailto:sales@trongenbio.com)

Website: [www.topscience.cn](http://www.topscience.cn)



## PRODUCT ADVANTAGES



Topscience adopts sterile production processes, eliminating the need for final sterilization steps, and ensuring product quality.



Ultra-low levels of impurities, bacterial endotoxins, residual elemental impurities, and microbial control.



Frost & Sullivan awarded Topscience Biotech Sterile Sodium Hyaluronate the honorary title of pioneer.

## PRODUCTION ADVANTAGES

Completely Closed  
Production Process



Tightly Controlled  
Production Environment



Advanced and Complete  
Production Equipments



Excellent Production  
Process Standards



## APPLICATIONS

Medicinal field:

Suitable for viscoelastic implantation, joint surgery, and other fillers, and can be used to prevent postoperative adhesions, etc.

Medical aesthetic field:

suitable for Class II and III medical devices and other fields.

## PACKAGE AND STORAGE

1. 100g/bottle, pharmaceutical aluminum bottle with outer pharmaceutical low-density polyethylene (LDPE) bag.
2. In an airtight container, protected from light and humidity, 2-8°C.